



# Agilent

**Errata Notice**

This document contains references to "Advanced Analytical" or "AATI." Please note that Advanced Analytical was purchased by Agilent in June 2018. For more information, contact Agilent via: [www.agilent.com/chem/contactus](http://www.agilent.com/chem/contactus)



## OLIGO PRO II



Automated Oligonucleotide Analyzer



# Accurately Assess Purity


## Ensure Quality Product

Oligonucleotide synthesis and purification is a complex, multi-step process. Final product quality is influenced by numerous factors that, while carefully monitored, cannot always be controlled. It is critical to both the manufacturer and end-user that oligonucleotide purity be known, and the purity be of sufficient quality for the end application, to ensure optimal performance and minimize ambiguous results.



*Purity analysis of oligonucleotides with the OLIGO PRO II Automated Oligonucleotide Analyzer. The electropherogram on the left shows a 57.99% pure oligonucleotide sample. The electropherogram on the right depicts a high purity oligonucleotide sample at 93.71%.*

Advanced Analytical's **OLIGO PRO II Automated Oligonucleotide Analyzer** provides high resolution, accurate assessments of oligonucleotide purity with unsurpassed sample throughput via parallel capillary gel electrophoresis. Capillary arrays with the option of 12, 24, or 96 capillaries are available and completely interchangeable, enabling seamless scalability. Streamlined operation, ready-to-use reagents, and dedicated software offer a complete solution for ensuring product quality.




*Broad Sizing Range - Patented Oligel separation matrix enables single nucleotide resolution through 60-mer lengths under appropriate sample load.*

# Award-Winning Technology

## Unique, Unmatched Capabilities

Advanced Analytical's parallel capillary gel electrophoresis-UV detection (CGE-UV) and separation technology dramatically improves laboratory throughput and efficiency while minimizing sample consumption and reducing operational costs.

CGE-UV is a standard method for assessing oligonucleotide purity, providing high separation resolution and direct detection of ssDNA and ssRNA without the sequence and size dependent variations encountered in mass spectrometry or fluorescence-based methods.



## Scalable

Quickly adjust throughput with interchangeable 12, 24, or 96 capillary arrays. Simple capillary array design is easy to handle.

## Direct Detection

Direct UV detection of oligonucleotides, no intercalating dye required.

## Sample Capacity

Load up to three 96-well plates for unattended oligonucleotide analysis. Plates can be loaded and prioritized while the system is running.

## Compact

Small instrument footprint saves bench space.



# Features and Benefits

- **Fast Separations**  
Complete 12, 24, or 96 separations in about 1 hour.
- **Improve Product Quality**  
Monitor synthesis performance to demonstrate purity of ssDNA or ssRNA.
- **Wide Sizing Range**  
Achieve n-1 resolution through 60-mers.
- **Throughput for All Labs**  
Variable throughput with the option of 12, 24, or 96 capillary arrays.
- **Direct Detection of Oligos**  
No dye labeling needed, uses UV detection.
- **Simplified Report Generation**  
Advanced peak integration software and data reporting capabilities.

## Specifications

**Array Sizes:** 12, 24, or 96 capillaries

**Detection:** Online, fixed wavelength UV absorbance at 254 nm

**Sample Injection:** Simultaneous electrokinetic injection from a 96-well microplate

**Sample Format:** Desalted to maximize injection efficiency; 1 – 5  $\mu$ M concentration range

**Sample Volume Required:** Typical volume 100  $\mu$ L/well (minimum of 20  $\mu$ L/well)

**Oligo Sizing Range:** Up to 60-mer lengths under appropriate sample load

**Software:** Proprietary **OLIGO PRO II** software for system control/data analysis

**Data Export Format:** CSV, PDF, JPEG, BMP, or PNG for individual samples or entire sample set

**Environmental Conditions:** Indoor use, normal laboratory environment; lab temperature 15 - 25°C

**Relative Humidity Range:** 15 – 60% (non-condensing)

**Electrical:** 100 – 200 VAC; 50 – 60 Hz (200 – 230 VAC; 50 – 60 Hz available); 15 A

**Instrument Dimensions:** 67.31 cm H x 48.26 cm W x 61 cm D (26.5 x 19 x 24 in)

**Instrument Weight:** 102 Kg (125 lbs)

[www.agilent.com](http://www.agilent.com)

For Research Use Only. Not for use in diagnostic procedures.

This information is subject to change without notice.

© Agilent Technologies, Inc. 2018  
Published in the USA, November 5, 2018  
5994-0421EN

OLIGO PRO II Automated Oligonucleotide Analyzer