

A Rapid Method for Bioburden Testing in Live Attenuated Influenza Vaccine Intermediates – From Concept to Implementation

Praful K. Bhusari, MedImmune Vaccines, Inc., Mountain View, CA 94043 and
Ann M. Steger, Advanced Analytical Technologies, Inc., Ames, IA 50010

Bioburden levels are routinely evaluated throughout the viral vaccine manufacturing process. A rapid method for monitoring bioburden in in-process influenza vaccine material has been developed through a collaborative effort between MedImmune Vaccines, Inc. (Mountain View, CA) and Advanced Analytical (Ames, IA). The live, attenuated influenza virus vaccine material is produced by inoculating fertile, specific pathogen-free (SPF) chicken eggs. Following incubation and propagation of the virus in the SPF eggs, the allantoic fluids are harvested, which have natural potential for harboring microbial contamination. One-liter sub-lots of allantoic fluids are collected and sequestered until bioburden test results are obtained, at which time the sub-lots with satisfactory bioburden results are pooled and further processed. The sub-lot quarantine and delay in bioburden results prevents a continuous process, thus warranting a rapid bioburden detection method. Initial feasibility and proof-of-concept were successfully completed on the *Micro PRO*[™] detection system using allantoic fluids spiked with *Staphylococcus aureus* and *E. coli*. Subsequently, a novel pre-treatment procedure was developed to significantly reduce background contribution from cellular debris and other allantoic fluid components. This collaboration has yielded a procedure for rapid screening of bioburden levels in in-process egg-based vaccine material.