

**Subject:** Removal of Conical Centrifuge Tubes from 900 Series Analysis Kits

**Effective Date:** July 2015

**Products Affected:** All Fragment Analyzer™ 900 Series Analysis Kits

Dear Customer,

We are contacting you to let you know that we have made a change to the kit components of the 900 Series Fragment Analyzer™ Analysis Kits. The 250 mL conical centrifuge tube that was previously supplied with the kit is no longer included.

The 50 mL or 250 mL size conical tubes used for solution preparation with the Fragment Analyzer™ can be ordered separately in bulk from several different vendors. Ordering information is provided below, and also in each kit User Manual.

Conical centrifuge tubes used for prepared Gel/Dye mixture and/or 1X Capillary Conditioning Solution:

- a. 250 mL (for 96-Capillary instruments or larger volumes): Corning #430776, available from Fisher #05-538-53 or VWR #21008-771
- b. 50 mL (for 12-Capillary instruments or 50 mL volumes): BD Falcon™ #352070, available from Fisher #14-432-22 or VWR #21008-940

All respective User Manuals have been updated to reflect this change in July of 2015.

Please contact AATI Technical Support with any questions; 515-296-6600 (current Ames site), 515-964-8500 (new Ankeny site) or [support@aati-us.com](mailto:support@aati-us.com).

Thank you,

AATI Customer Service Team