

Subject: Addition of FS-SLR905 35-400 bp Range DNA Ladder to DNF-900 Kits

Effective Date: October 2015

Products Affected:

1. DNF-900 DNF-900 dsDNA Reagent Kit, 35 bp - 500 bp
 - a. FS-SLR905-0001 35-400 bp Range DNA Ladder, 1.0 mL (**Added**)
 - b. FS-SLR900-0001 75-400 bp Range DNA Ladder, 1.0 mL (**Removed/Replaced**)

A. Kit Components

1. Replaced 75-400 bp Range DNA Ladder (# FS-SLR900) with 35-400 bp Range DNA Ladder (# FS-SLR905) in the DNF-900 kit assemblies. The new ladder now includes a 35 bp DNA fragment, and is provided at a total ladder concentration of 2.0 ng/ μ L.
2. The 75-400 bp Range DNA Ladder (# FS-SLR900) is being discontinued due to supplier changes, effective October 2015.
3. These changes do not affect the performance of the DNF-900 kits.

All respective User Manuals have been updated to reflect this change in October of 2015.

Please contact AATI Technical Support with any questions; 515-964-8500 or support@aati-us.com.

Thank you,

AATI Customer Service Team