

**Subject:** Addition of FS-SLR905 35-400 bp Range DNA Ladder to DNF-912 Kits

**Effective Date:** June 15, 2016

**Products Affected:**

1. DNF-912 dsDNA 912 Reagent Kit, 35 bp - 500 bp
  - a. FS-SLR905-0001 35-400 bp Range DNA Ladder, 1.0 mL (**Added**)
  - b. FS-SLR900-0001 75-400 bp Range DNA Ladder, 1.0 mL (**Removed/Replaced**)

A. Kit Components

1. Replaced 75-400 bp Range DNA Ladder (# FS-SLR900) with 35-400 bp Range DNA Ladder (# FS-SLR905) in the DNF-912 kit assemblies. The new ladder now includes a 35 bp DNA fragment, and is provided at a total ladder concentration of 2.0 ng/ $\mu$ L.
2. The 75-400 bp Range DNA Ladder (# FS-SLR900) is being discontinued due to supplier changes.
3. These changes do not affect the performance of the DNF-912 kits.

All respective User Manuals will be updated to reflect this change by June 30, 2016.

Please contact AATI Technical Support with any questions; 515-964-8500 or [support@aati-us.com](mailto:support@aati-us.com).

Thank you,

AATI Customer Service Team