

Subject: Extension of Shelf Life of DNF-375-0120 High Sensitivity Genomic DNA Diluent Marker Solution

Effective Date: October 5, 2016

Products Affected:

1. DNF-467-0500 Genomic DNA 50Kb Analysis Kit, 500 Samples
 - a. DNF-375-0120 High Sensitivity Genomic DNA Diluent Marker Solution, 120 mL (**Shelf Life Extended from 6 months to 12 months after manufacture**)

A. Extension of Product Shelf Life

1. The shelf life of DNF-375-0120, High Sensitivity Genomic DNA Diluent Marker Solution, 120 mL, has been extended from 6 months to 12 months from the date of manufacturing.
2. This change was made following internal evaluations, which demonstrated equivalent performance of the product when stored under the recommended conditions of 2-8°C temperatures, for varying periods of time up to and exceeding 12 months after the date of manufacturing.

All respective packaging and documentation going forward will be updated to reflect this change by October 5, 2016.

Please contact AATI Technical Support with any questions; 515-964-8500 or support@aati-us.com.

Thank you,

AATI Customer Service Team